

**General Advice:** These instructions are general advice only and do not take into account your specific circumstance. These instructions are to be used in conjunction with the information provided by your health professional and any instructions provided by the equipment manufacturer. Do not exceed the Maximum User Weight specified by the manufacturer. If you are unsure on how to use or fit a piece of equipment, please contact your Health Professional or Country Care Store.

## Bed & Chair Raisers Use

- Bed and chair raisers are used to raise the height of a bed or chair, making it easier for people with mobility issues to get in and out of the bed or chair. They are typically used by people with conditions such as arthritis, back pain, or other conditions that make it difficult to bend or stand from low heights
- Bed and chair raisers come in different height sizes and materials, such as wood, plastic, and metal and are applied under the legs of a chair or bed. Raisers can have different size indentations to accommodate the various dimensions of the chair or bed legs. Before purchase ensure the chair or bed leg dimensions will fit in the raiser indentations

## To Install

1. Place the raisers under the legs of the bed or chair, making sure that they are level and stable, you may need someone to assist you
- **Note:** Make sure that the weight capacity of the raisers is compatible with the weight of your furniture
  - Beds occasionally have a number of legs. A raiser will need to be put under each leg to ensure stability
  - If moving the furniture at any time ensure the legs of the chair or bed are still securely in the raisers before use

## PRECAUTIONS & SAFETY

- Some beds have more than 4 legs. You will need to have one raiser under each leg of the bed or chair to ensure its stability. Thus a bed with 6 legs will require 6 raisers
- If you move your bed or chair to a new position it is important to ensure that all chair or bed legs/castors are held securely in the raisers. Small movements through regular use may also alter the alignment of the raisers

## MAINTENANCE & CLEANING

- Please check the product label / user manual for specific cleaning instructions
- If the equipment appears damaged, **DO NOT USE**. It may require replacing



*Chair / Bed Raisers also available in cone shape*