



# Swivel Bather & Bath Chairs

## INSTALLATION & SAFE USE GUIDE

These instructions are to be used in conjunction with the information that your Occupational Therapist has already discussed with you.

### What is a Swivel Bather Used For?

- A swivel bather enables anyone who has decreased balance, physical strength or mobility difficulties to bathe safely and independently while sitting over the bath.
- A swivel bather has a seat, backrest and arms. The seat is secured to a base frame which fits securely onto the bath. It also incorporates a locking system to enable safe positioning of the seat.

### TO INSTALL:

#### Fitting the base frame over the bath:

- Rest the base of the frame securely on the bath.
- Turn the positioning screws on both sides until the black rubber ends are against the bath. This will prevent the base from moving or slipping.

#### Fitting the seat:

- Place the seat on the base frame. The socket under the seat locates over the swivel pin of the base frame. Ensure that the seat rotates freely.

### How to Use a Swivel Bather

- Lock the seat by pushing the knob on the locking lever gently down until the lock locates in the slots in the base frame.
- Swivel the seat to face the side of the bath and lock into position. The user can now be seated.
- Unlock the seat by lifting the locking lever and transfer legs over the side of the bath.
- Lock the seat when over the bath for extra security while washing or showering.



### PRECAUTIONS & SAFETY

- Ensure that the seat is locked into position before sitting or standing.
- Ensure that the items required for showering are in easy reach when seated.

**TO BE USED AS PER THE MANUFACTURERS WEIGHT RECOMMENDATIONS**

### MAINTENANCE & CLEANING

- Ensure that you regularly clean the chair with a mild disinfectant/cleanser.
- Should the equipment appear damaged – **DO NOT USE**. It may require replacing.